

International Conference on Public Health, Epidemiology & Infectious Diseases (ICPEI 2026)

Conference Proceedings

June 25-26, 2026 - Taipei, Taiwan

Serum Determinations in SARS COV2 Infection

Dr. arch.eng. Bârdan Marius Vasile

Research and Development Center “Dr. Bârdan Marius Vasile” TIC 50890649, Romania

Abstract

Purpose

This presentation outlines the technical specifications, procedural guidelines, and clinical performance data for the VIDAS® SARS-CoV-2 IgG II (9COG) assay (pp. 1-2). The test serves as an automated diagnostic aid to evaluate whether an individual has developed a specific humoral immune response following exposure to the SARS-CoV-2 virus or vaccination (pp. 2, 4).

Methodology

The assay utilizes the Enzyme Linked Fluorescent Assay (ELFA) technique, combining a two-stage sandwich enzyme immunoassay with automated fluorescence detection via a specialized Solid Phase Receptacle (SPR) device (pp. 2, 5). The inner wall of the SPR is coated with recombinant receptor-binding domain (RBD) of the viral Spike protein to capture specific IgG antibodies from human serum or plasma (pp. 2, 7). Alkaline phosphatase-labeled anti-human IgG antibodies and a 4-methylumbelliferyl phosphate substrate are then introduced to generate a measurable fluorescent signal (pp. 6, 11).

Performance and Results

The automated analysis delivers semi-quantitative index results within approximately 27 minutes (pp. 2, 24). The test is highly precise and exhibits a wide linear range between 0.40 and 17.00 index values, with a determined cut-off value of 1.00 index equivalent to 20.33 BAU/ml relative to the first WHO international standard (pp. 27, 29, 31). Clinical evaluations demonstrate an overall internal specificity of 99.9% (p. 39). Clinical sensitivity is strongly correlated with time elapsed post-infection, rising from 45.3% within the first 7 days following a positive PCR result to 96.6% after 16 or more days (p. 40). Cross-reactivity studies against a panel of 261 potentially interfering infectious conditions verified high analytical specificity with minimal cross-reactivity (pp. 35-36).

Conclusion

The VIDAS® SARS-CoV-2 IgG II assay provides a highly efficient, reliable, and standardized solution for clinical laboratories to quantify anti-RBD Spike protein IgG antibodies, facilitating accurate long-term humoral immunity assessment (pp. 2, 4).

Keywords: SARS-CoV-2, Immunoglobulin G (IgG), Enzyme Linked Fluorescent Assay (ELFA), Receptor Binding Domain (RBD), Spike Protein, VIDAS Automation, Humoral Immunity, Clinical Sensitivity, Analytical Specificity